


MODEL DH29-46SD

45 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT		WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	48 IN	25.75 IN	62 IN		121 CM	65 CM	157 CM
EXTERIOR DIMENSIONS	52 IN	34 IN	82 IN		133 CM	87 CM	209 CM
CRATED DIMENSIONS	53 IN	43 IN	86 IN		134 CM	109 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 475 LBS.

SHIPPING WEIGHT: 525 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)

DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)

INSIDE DELIVERY (WHITE GLOVE DELIVERY)