


MODEL DH4-80PT

80 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	75 IN	30 IN	59 IN	190 CM	76 CM	149 CM
EXTERIOR DIMENSIONS	78 IN	37.5 IN	81.5 IN	198 CM	95 CM	207 CM
CRATED DIMENSIONS	85 IN	41 IN	86 IN	215 CM	104 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 800 LBS.
SHIPPING WEIGHT: 850 LBS.

SHIPPING OPTIONS

- DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)
- DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)
- INSIDE DELIVERY (WHITE GLOVE DELIVERY)