


MODEL DHK4-14GD

14 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	22.75 IN	21 IN	53.25 IN	58 CM	53 CM	135 CM
EXTERIOR DIMENSIONS	25.75 IN	26.5 IN	78.5 IN	65 CM	67 CM	199 CM
CRATED DIMENSIONS	32 IN	32 IN	80 IN	81 CM	81 CM	203 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 190 LBS.

SHIPPING WEIGHT: 240 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)

DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)

INSIDE DELIVERY (WHITE GLOVE DELIVERY)