


MODEL DHK4-49GD

46.62 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


| | WIDTH | DEPTH | HEIGHT | WIDTH | DEPTH | HEIGHT |
|---------------------|----------|----------|----------|--------|--------|--------|
| INTERIOR DIMENSIONS | 52.75 IN | 26.25 IN | 59.5 IN | 134 CM | 67 CM | 152 CM |
| EXTERIOR DIMENSIONS | 57.25 IN | 34.75 IN | 78.75 IN | 145 CM | 88 CM | 200 CM |
| CRATED DIMENSIONS | 85 IN | 41 IN | 86 IN | 215 CM | 104 CM | 218 CM |

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW

AUTHORIZED DELAER OF **LIEBHERR**

WWW.SO-LOW.COM



NET WEIGHT : 750 LBS.

SHIPPING WEIGHT: 800 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)

DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)

INSIDE DELIVERY (WHITE GLOVE DELIVERY)