

# MODEL DHP4-74SD


72 CUBIC FOOT LAB & PHARMACY REFRIGERATOR

	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	72 IN	27 IN	55 IN	183 CM	69 CM	140 CM
EXTERIOR DIMENSIONS	75 IN	32 IN	81 IN	190 CM	81 CM	206 CM
CRATED DIMENSIONS	85 IN	41 IN	86 IN	215 CM	104 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

## SIDE VIEW

## FRONT VIEW



**SO-LOW**

[WWW.SO-LOW.COM](http://WWW.SO-LOW.COM)



**NET WEIGHT : 690 LBS.**  
**SHIPPING WEIGHT: 720 LBS.**

### SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )  
DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )  
INSIDE DELIVERY ( WHITE GLOVE DELIVERY )