


# MODEL LRT30W1HC

## 29.1 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	24.5 IN	31.75 IN	59 IN	62 CM	80 CM	150 CM
EXTERIOR DIMENSIONS	31.25 IN	38.75 IN	85.25 IN	79 CM	98 CM	216 CM
CRATED DIMENSIONS	46 IN	43 IN	86 IN	116 CM	109 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

### SIDE VIEW



### FRONT VIEW



# SO-LOW

AUTHORIZED DELAER OF **LIEBHERR**

[WWW.SO-LOW.COM](http://WWW.SO-LOW.COM)



**NET WEIGHT : 350 LBS.**

**SHIPPING WEIGHT: 400 LBS.**

### SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )

DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )

INSIDE DELIVERY ( WHITE GLOVE DELIVERY )