


MODEL U85-25X

25 CUBIC FOOT ULTRA LOW TEMPERATURE FREEZER


	WIDTH	DEPTH	HEIGHT		WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	34 IN	25 IN	51.25 IN		86 CM	63 CM	130 CM
EXTERIOR DIMENSIONS	54.75 IN	37 IN	79.5 IN		139 CM	93 CM	201 CM
CRATED DIMENSIONS	68 IN	45 IN	86 IN		172 CM	114 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



*TO PASS THROUGH DOORWAYS, THE DOOR CAN OPEN A FULL 180°, TO REDUCE THE UNIT EXTERIOR DEPTH TO 33"

SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 900 LBS.

SHIPPING WEIGHT: 1125 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)

DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)

INSIDE DELIVERY (WHITE GLOVE DELIVERY)