

MODEL VSC40-22


22 CUBIC FOOT ULTRA LOW TEMPERATURE FREEZER

| | WIDTH | DEPTH | HEIGHT | WIDTH | DEPTH | HEIGHT |
|---------------------|--------|-------|---------|--------|--------|--------|
| INTERIOR DIMENSIONS | 67 IN | 19 IN | 30 IN | 170 CM | 48 CM | 76 CM |
| EXTERIOR DIMENSIONS | 98 IN | 36 IN | 43.5 IN | 248 CM | 91 CM | 110 CM |
| CRATED DIMENSIONS | 101 IN | 43 IN | 53 IN | 256 CM | 109 CM | 134 CM |

NOTE THAT THE FREEZER LID CAN OPEN UP TO APPROXIMATELY 75° FROM HORIZONTAL POSITION

SIDE VIEW

FRONT VIEW



SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 850 LBS.
SHIPPING WEIGHT: 1075 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)
 INSIDE / WHITE GLOVE DELIVERY (REQUIRES A SEMI-TRUCK RECEIVING DOCK)