


# MODEL VSU40-25

## 25 CUBIC FOOT ULTRA LOW TEMPERATURE FREEZER


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	34 IN	25 IN	51.25 IN	86 CM	63 CM	130 CM
EXTERIOR DIMENSIONS	50 IN	37 IN	79.5 IN	127 CM	93 CM	201 CM
CRATED DIMENSIONS	53 IN	43 IN	86 IN	134 CM	109 CM	218 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

**SIDE VIEW**



**FRONT VIEW**



\*TO PASS THROUGH DOORWAYS, THE DOOR CAN OPEN A FULL 180°, TO REDUCE THE UNIT EXTERIOR DEPTH TO 33"

**SO-LOW**

[WWW.SO-LOW.COM](http://WWW.SO-LOW.COM)



**NET WEIGHT : 900 LBS.**  
**SHIPPING WEIGHT: 1125 LBS.**

**SHIPPING OPTIONS**

- DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )
- DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )
- INSIDE DELIVERY ( WHITE GLOVE DELIVERY )