


MODEL VSU80-30

30 CUBIC FOOT ULTRA LOW TEMPERATURE FREEZER


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	42 IN	25 IN	51.25 IN	106 CM	63 CM	130 CM
EXTERIOR DIMENSIONS	60.5 IN	37 IN	79.5 IN	153 CM	93 CM	201 CM
CRATED DIMENSIONS	68 IN	45 IN	86 IN	172 CM	114 CM	218 CM

NOTE THAT SHELVES ARE FIXED AND CANNOT BE REMOVED OR ADJUSTED UNLESS SPECIAL ORDERED AS ADJUSTABLE.

SIDE VIEW



FRONT VIEW



*TO PASS THROUGH DOORWAYS, THE DOOR CAN OPEN A FULL 180°, TO REDUCE THE UNIT EXTERIOR DEPTH TO 33"

SO-LOW

WWW.SO-LOW.COM



NET WEIGHT : 975 LBS.

SHIPPING WEIGHT: 1325 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)

DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)

INSIDE DELIVERY (WHITE GLOVE DELIVERY)