

## MODEL DHK30-24MDP

### MANUAL DEFROST FREEZER -15°C TO -25°C

<b>Capacity [ Cu. Ft. ]</b>	<b>23.5 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>668 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-15°C to -30°C</b>
<b>Temperature Range [ °F ]</b>	<b>5°F to -22°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	34.75" x 37" x 76"
Interior [ Inches ]	27.5" x 28" x 51"
Exterior [ Centimeters ]	88 x 94 x 193
Interior [ Centimeters ]	70 x 71 x 129

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R290
Exterior	Powder Coated Steel
Interior	Aluminum Liner
Insulation	Polyurethane Foam
Lockable Door	Keyed
Porthole	1" Diameter
Casters	4, with 2 Leveling Feet
Storage	4 Adjustable Shelves
Shipping Weight	215 lbs.


#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure
Data Logging	Yes, Adjustable Intervals
Data Download	Yes, via USB, PDF Format
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts
Backup Battery	Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor