

MODEL DHK4-14GD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 14 Cu. Ft.
Capacity [Liters] 395 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 35°F to 46°F

Dimensions Width x Depth x Height
Exterior [Inches] 25.75" x 26.5" x 78.5"
Interior [Inches] 23" x 21" x 53.25"
Exterior [Centimeters] 65 x 67 x 199
Interior [Centimeters] 58 x 53 x 135

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Automatic
Refrigerant R600a
Exterior Steel
Interior Polystyrene (HIPS)
Insulation US EPA and SNAP approved
Lockable Door Yes, Keyed
Porthole Yes, 1" Diameter
Casters 4, with Leveling Feet
Storage 6 Adjustable Shelves, 1 Basket
Shipping Weight 250 lbs.


Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 3.77A
Power Consumption 3.71 kW/d
NEMA Configuration NEMA 5-15
Certification(s) UL Listed

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp, Door Ajar
Power Failure, Low Battery
Sensor Failure, USB Failure
Data Logging Yes, Adjustable Intervals
Data Download Yes, via USB, PDF Format
Min / Max Temperature Yes, Display and Reset
Alarm Relay Yes, Dry Contacts
Backup Battery Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor