

SHELDON MODEL SRI3P

REFRIGERATED INCUBATOR

Capacity [Cu. Ft.]	3.5 Cu. Ft.
Capacity [Liters]	99 Liters
Temperature Range [°C]	15°C to 40°C (At 25°C ambient)
Temperature Range [°F]	59°F to 104°F (At 77°C ambient)

Dimensions	Width x Depth x Height
Exterior [Inches]	24.25" x 21.5" x 34"
Interior [Inches]	18.75" x 16.75" x 26.25"
Exterior [Centimeters]	61 x 54 x 85
Interior [Centimeters]	48 x 42 x 67

Air Circulation Forced Air

Exterior	Steel
Interior	Aluminum
Exterior Door	Steel
Interior Outlet	1
Porthole	1, Rear
Storage	Shelves (2 Standard, 3 Max)
Shipping Weight	125 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 4.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed CSA EN CE

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp, Door Ajar
Temperature Uniformity	+/-0.5°C (Midrange)
Temperature Stability	+/- 0.1°C (Midrange)

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Years Parts & Labor.