

## MODEL U85-22

### UPRIGHT ULTRA LOW FREEZER -40°C TO -85°C

**Capacity [ Cu. Ft. ]** 22 Cu. Ft.  
**Capacity [ Liters ]** 626 Liters  
**Temperature Range [ °C ]** -40°C to -85°C  
**Temperature Range [ °F ]** -40°F to -120°F

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

**Dimensions**  
 Width x Depth x Height  
 Exterior [ Inches ] 43" x 37" x 79.5"  
 Interior [ Inches ] 28.5" x 26" x 51.25"  
 Exterior [ Centimeters ] 109 x 93 x 201  
 Interior [ Centimeters ] 72 x 66 x 130

**Refrigeration System**  
 Cascade-Style System  
 Compressors 2, Hermetic Compressors  
 Air Circulation Gravity Flow  
 Defrost Method Manual  
 Refrigerant HCFC & CFC Free  
 Exterior Powder Coated Steel  
 Interior Zinc Coated Galvanized Steel  
 Insulation Polyurethane (Ecomate) Foam  
 Inner Doors 5, Stainless Steel  
 Lockable Door Lockable with Padlock  
 Porthole Optional 1/2"  
 Casters 4, Swivel with Locks  
 Storage 4 Adjustable Shelves, with 1 Base Shelf

Shipping Weight 975 Lbs.



**Electrical**  
 Standard Voltage 115 Volts, 60 Hertz, 1 Phase  
 Amps / Run Amps / Plug 20A | 15A | NEMA 5-20  
 ( Optional Voltages ) 208 Volts, 230 Volts  
 Certification(s) ETL Listed

**Controller**  
 Display Type Microprocessor, Digital  
 Alarms High/Low Temp, Power Failure, Low Battery  
 Alarm Relay Yes, Dry Contacts  
 Automatic Restart Yes, Automatically Restarts After Power Failure  
 Backup Battery Yes, Only Powers Alarm



**Warranty** 1 Year Parts & Labor

Specifications subject to change without notice